



# 008, Wing-3, Devi Krupa,  
Pant Nagar, Ghatkopar (E)  
Mumbai - 400 075

E-mail : [contact@crystaldatainc.com](mailto:contact@crystaldatainc.com)  
[crystaldatainc@gmail.com](mailto:crystaldatainc@gmail.com)  
[www.crystaldatainc.com](http://www.crystaldatainc.com)



NSIC-CRISIL SE3C

# Crystal Data Inc.

CONSULTING, DEVELOPMENT, SUPPORT

Phone : 022-21022438 / 21021849 / 65883988  
Mobile : 9820373936 / 9920548030

LAB NO. : 5 DATE : 12-Aug-2011  
PATIENT NAME : MR. KETAN CHAVAN SEX : Male  
REF. BY DR. : DR. PATIL M.B.B.S. AGE : 29 Years  
SAMPLE COLL. AT : CRYSTAL LAB

## SLIDE TEST FOR ANTI - STREPTOLYSIN O

Result Positive  
Titre 1:120  
Method Latex Agglutination  
Kit Used RHELAX - ASO TULIP  
Control Used Positive & Negative

Comments : Rhelax ASO slide test for detection of antibodies to streptolysin O is based on the principle of agglutination. The test specimen (serum) is mixed with Rhelax ASO latex reagent and allowed to react. If antibodies to streptolysin O are present in concentrations more than 200 IU/ml but less than 400 IU/ml then a visible agglutination is observed. If antibodies to streptolysin O are not present or are in concentration less than 200 IU/ml then no agglutination will be observed.

Agglutination is a positive test result and indicates the presence of detectable levels of Anti-streptolysin O in the test specimen. No agglutination is a negative test result and indicates the absence of detectable levels of Anti-streptolysin O in the test specimen.

----- End of Report. -----

Chandan Vartak  
D.M.L.T.

Highlighted Result Values Indicate Abnormal

Report Printed By My Lab [www.crystaldatainc.com](http://www.crystaldatainc.com)

Dr. Pankaj Shah  
M.D. M.B.B.S.